ASTAR. physiotechnology

PhysioGo 701C o -













It is an unit with 3 independent treatment channels for electrotherapy, ultrasound therapy, combined therapy and laser therapy. In addition,

- 379 preset treatment programs
- 38 treatment sequences for electrotherapy
- 350 user-defined programs
- legible encyclopedia format allow the therapist to utilize its full potential.
- The battery allowing for work up to 8 hours without recharging

Product Background:

- made in Poland (Europe)
- more than 20 years of experience
- with CE certificate















PhysioGo 701C Specifications



7" color touchscreen



Number of preset treatment sequences



Operation by the keyboard



Number of sequences to be set by the users



Three completely independent treatment channels



Electrodiagnostic of the nervous-muscle system (I/t)



Encyclopedia



Electrode test



Device weight, dedicated bag



Number of currents generated by the device



Battery in standard

Number of preset

treatment programs



Number of multi-frequency (1; 3,5 MHz) ultrasound heads which the device supports




Number of laser applicators which the device supports



Number of programs to be set by the users







Both 4cm² and 1 cm² ultrasound head offer multi-frequency which has both 1MHz for deep area and 3.5MHz for superficial area.



Lithium-ion Battery allows operation up to 8 hrs. Cable free mode increases mobility practice



The feet of the device lets you change the height for easy operation.



www.tuv.com ID 0000034043



QUALITY OF THE YEAR 2010 – GOLD



Download PDF





EN ISO 13485:2016 Fulfilled the ISO quality management system

www.newtrendmedic.com/physiogo



i-NOVO AWARDS 2015 PhysioGo





























